

LASAK BioniQ QR Exocad library structure

LASAK_BioniQ_(Ti Base).zip

Scanbody QR, indexed – long
Ref. No. 2815.00

Scanbody Screw-On – long
Ref. No. 2835.00


Scanbody, CEREC
indexed – Bluecam/L
Ref. No. 2821.00

Screw_QR

Dummy base shown
in the Exocad library


The real base of the
BioniQ component

bioniq_qr_interface




FITS
to

BioniQ QR Ti base
L0.8




Ref. No.
2159.00

bioniq_qr_L2.0_interface




FITS
to

BioniQ QR Ti base
L2.0




Ref. No.
2159.20

bioniq_qr_ni_interface




FITS
to

BioniQ QR/NI Ti base
L0.8




Ref. No.
2186.00

bioniq_qr_ni_L2.0_interface




FITS
to

BioniQ QR/NI Ti base
L2.0




Ref. No.
2186.20

bioniq_qr_b_interface




FITS
to

BioniQ QR/B Ti base
L0.4



Ref. No.
2158.00

bioniq_qr_b_L2.0_interface




FITS
to

BioniQ QR/B Ti base
L2.0




Ref. No.
2158.20

bioniq_so46_ni




FITS
to

BioniQ Screw-On




Ref. No.
2148.01
2148.02
2148.03
2148.04
2149.03
2149.04
2149.05
2150.04
2150.05

bioniq_cerec_interface



FITS
to

BioniQ QR, Ti base CEREC




Ref. No.
2183.00

Free library

LASAK_BioniQ_(Ti Base).zip

Ref. No.
2858.00




S3.5
S4.0
T4.0
S5.0
T5.0

LASAK BioniQ QN Exocad library structure


LASAK_BioniQ_(Ti base).zip

Scanbody QN, indexed – long
Ref. No. 2840.00




Dummy base shown
in the Exocad library


bioniq_qn_interface




bioniq_qn_L2.0_interface




bioniq_qn_ni_interface



bioniq_qn_ni_L2.0_interface




BioniQ QN Ti base
L0.8




Ref. No.
2181.00

BioniQ QN Ti base
L2.0




Ref. No.
2181.20

BioniQ QN/NI Ti base
L0.8




Ref. No.
2189.00

BioniQ QN/NI Ti base
L2.0



Ref. No.
2189.20

bioniq_so46_ni




BioniQ Screw-On




Ref. No.
2177.01
2177.02
2177.03
2177.04

bioniq_cerec_interface




BioniQ QN, Ti base CEREC



Ref. No.
2198.00


Scanbody Screw-On – long
Ref. No. 2835.00



Scanbody CEREC,
indexed – Bluecam/S
Ref. No. 2864.00



Screw_QN




The real base of the
BioniQ component

Free library

LASAK_BioniQ_(Ti base).zip

Ref. No.
2859.00



S2.9

